



# Safety Waste Sets Overview

Sets for the safe disposal of liquid laboratory waste.

- ✓ Made from durable and robust PE-HD material.
- ✓ Sets with a practical lid or solid ball valve funnel.
- ✓ All canisters with UN-Y approval for safe transport and storage.
- ✓ Handy canister volume from 2.5 to 10 litres - ideal for laboratory applications.
- ✓ All sets can be used immediately without additional accessories.

All set components can also be reordered individually:



Component	Funnel		Canister			
	Ball valve	Hinged lid	Safety canister	Safety canister	Safety canister	Safety canister
Type	Ball valve	Hinged lid	Safety canister	Safety canister	Safety canister	Safety canister
Model	<b>ARNOLD</b>	<b>Marco</b>	-	-	-	-
Feature	with lance (220 mm)	without lance	-	-	-	-
Content	-	-	2.5 litre	5 litre	10 litre	10 litre
Thread	S55	GL45	GL45	GL45	GL45	S55
Set	Part number					
A	307 316	•	•			
B	307 317			•		
C	307 318		•		•	
D	307 450	•				•